Endostar Endosizer

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Endostar Endosizer

Instructions for cleaning, disinfection and sterilisation

The supplied instrument is non-sterile. Before the first use and after each subsequent use, the instrument must be cleaned, disinfected and sterilised in accordance with the following instructions.

1. Recommendations

- The instrument can be disinfected using mild spray or foam disinfectants.
- Disinfectants must be suitable for disinfecting instruments made of plastic, free from aldehydes, phenol, anti-corrosives (especially di- or triethanolamine), oxidants (hydrogen peroxide, sodium hypochlorite > 5%), solvents.
- Please use personal protective equipment for your own safety, e.g. gloves, safety glasses, a face mask and a face shield.

2. Warnings

- When using chemical soaking solutions, always follow the liquid manufacturer's instructions.
- Excessive fluid levels and excessive soaking times (incorrect in accordance with the fluid manufacturer's instructions) may damage the instrument.
- Never expose the instrument to temperatures above 142°C (288°F). Make sure that the instrument is not placed in the autoclave too close to the walls or the bottom of the steam steriliser (the risk of excessive temperature and deformation).

INSTRUCTIONS FOR CLEANING, DISINFECTION AND STERILISATION OF THE INSTRUMENT

1. Disinfection

Soak the instrument or apply a liquid product with a spray tip immediately after use. Avoid soaking for long periods of time.

2. Washing

Rinse the instrument thoroughly using running tap water. The water should have room temperature.

3. Cleaning

Gently clean all surfaces of the instrument by immersing them in water with detergent. Use only soft brushes made of nylon, polypropylene or acrylic. After cleaning, rinse the instrument several times in clean water, preferably recently boiled, until the foam disappears. Use distilled water for the final rinse. Surfaces of the instrument should be visibly clean and free of stains.

4. Washing

Rinse the instrument thoroughly using demineralised water. The water should have room temperature.

5. Drying

Dry the instrument with a non-fibre cloth. Dry the instrument until there are no visible traces of moisture. The instrument should be thoroughly dried before packing.

6. Packaging

Place the instrument in a paper and foil bag that is suitable for steam sterilisation.

7. Sterilisation

Place the bag in the steriliser following the steriliser manufacturer's recommendations. If several instruments are sterilised in one autoclave cycle, ensure that the maximum load on the steriliser is not exceeded. Sterilisation in a steam autoclave using a pre-vacuum cycle is recommended. The instrument should be sterilised at 134°C for at least 4 minutes and at an overpressure of 2.1 bar or at 121°C for at least 20 minutes.

- The effectiveness and acceptance criteria of the procedure should be monitored.
- Please keep records regarding all processes and determine the shelf-life date following the packaging manufacturer's instructions.

8. Storage

Store the sterile-packaged instrument in a well-ventilated room, protected from moisture, dust, sunlight and heat. Before use, please visually assess that the integrity of the kit has not been compromised.