



Indications & Materials*

- Metal/Metal based
- Glass Ceramics (e.g. e.max, VITA Mark-II, IPS Empress 2)
- Oxide Ceramics (e.g. Lava™ Plus, Brux Zir, Procera)
- Resin Composites
- Resin Nano Ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)

*Only primary cement recommendations are shown. See respective instructions for use for full list of indications.

Cement Selection Guide

3M Science.
Applied to Life.™

3M Oral Care, 2510 Conway Avenue, St. Paul, MN 55144, USA



Veneers

- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ Veneer
- 3M™ RelyX™ Veneer
- 3M™ RelyX™ Ultimate

Inlays/Onlays

- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ Ultimate

Bridges

- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ Ultimate

Crowns

- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate

Endodontic Posts

- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Luting 2

Maryland Bridges

- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ U200

On Implant Abutments

- 3M™ RelyX™ U200
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ Temp NE
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200
- 3M™ RelyX™ Luting 2
- 3M™ RelyX™ Temp NE
- 3M™ RelyX™ Ultimate
- 3M™ RelyX™ U200

+



Features & Benefits

Cement Selection Guide

3M Science.
Applied to Life.™

3M, ESPE, Lava and ReQ-It are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies. © 2016. All rights reserved. 70-2013-0222-1

3M™ RelyX™ Temp NE/E
The Temporary Cement



- Available in both eugenol and non-eugenol formulation
- Easy removal of provisional
- Quick, easy excess removal
- Reliable, temporary adhesion

3M™ RelyX™ Ultimate
The CAD/CAM and Glass Ceramic Cement



- World-class esthetics with natural fluorescence
- Total-etch, self-etch or selective etch with just one adhesive
- Fast and easy procedure with only two components and few working steps
- Industry-leading bond strength



The Veneer Cement
3M™ RelyX™ Veneer

- Customer-preferred shading system
- Extremely color stable
- Matching try-in pastes
- Light-cure formula
- Easy handling and clean-up

+

The PEM and Pediatric Crown Cement
3M™ RelyX™ Luting 2



- Tack light-cure feature for "clean-up on demand"
- Easy to use automatic delivery
- Sustained fluoride release
- Virtually no post-operative sensitivity

- Saves time by eliminating etching, priming and bonding steps
- Virtually no post-operative sensitivity
- Strong and moisture-tolerant
- Endo tips make cementing posts easier and more reliable
- The world's most clinically proven self-adhesive resin cement



The All-Rounder

3M™ RelyX™ U200